

## **CGM Comparison Tool**

DEVICE	Dexcom G5	Dexcom G4	Guardian Connect	Guardian Link 2	MiniLink
Pairs with	<ul> <li>Specified mobile Devices         www.dexcom.         com/ous-         compatibility-         page</li> <li>Dexcom receiver         (not included in subsidy)</li> </ul>	<ul> <li>Dexcom receiver (not included in subsidy)</li> <li>Animas Vibe Pump</li> </ul>	Specified mobile devices www.medtronic-diabetes.com. au/guardian-connect	Medtronic 640G pump	Medtronic Veo pump
Can give insulin with CGM glucose reading without a finger prick	YES	NO	NO	NO	NO
Allows for up to 5 "followers"	YES	NO	YES	NO	NO
Low glucose suspend function	NO	NO	NO	YES	YES
Predictive low glucose suspend function	NO	NO	NO	YES	NO
Transmitter warranty life	3 months	6 months	12 months	12 months	12 months
Sensor life	7 days	7 days	6 days	6 days	6 days



## **CGM Comparison Tool**

DEVICE	Dexcom G5	Dexcom G4	Guardian Connect	Guardian Link 2	MiniLink
Device range	6m	6m	6m	2m	6m
Calibrations required	Once every 12 hours after initial set up	Once every 12 hours after initial set up	Once every 12 hours after initial set up	Best results: 3 - 4 times/day	Best results: 3 - 4 times/day
Accuracy is affected by Panadol	YES	YES	NO	NO	NO
Can over tape be used on the transmitter?	NO	NO	YES	YES	YES
What is the cost to me if using the subsidy?	<ul> <li>Cost of compatible devices if you do not have this already</li> <li>Cost of receiver if you would like to use this (not essential)</li> </ul>	Cost of the receiver if you would like to use this (not essential if you have an Animas pump)	Oost of compatible devices if you do not have this already	N/A	N/A

## For more information:

Dexcom G4: http://amsldiabetes.com.au/item/dexcom-g4-platinum/

Dexcom G5: http://amsldiabetes.com.au/item/dexcom-g5-mobile/

Medtronic Guardian Connect: https://www.medtronic-diabetes.com.au/guardian-connect

Medtronic Guardian Link 2 and MiniLink: https://www.medtronic-diabetes.com.au/pump-therapy/continuous-glucose-monitoring